



merits

SCOOTER COMPARISON

								
	Yoga	Roadster	Indigo	Urban	Cross Country*	Silverado Extreme	Vision Sport Power Chair**	Regal EZ Power Chair
Backrest Width	405mm	460mm	460mm	510mm	560mm	510mm	460mm	460mm
Seat Depth	350mm	460mm	410mm	460mm	500mm	460mm	460mm	400mm
Ground to Seat Height	595mm	650-700mm	640-690mm	680-730mm	560-640mm	600-680mm	560-750mm	520-560mm
Backrest Height	315mm	510mm	460mm	510mm	560mm	630mm	510mm	480mm
Overall Width	435mm	640mm	630mm	660mm	710mm	720mm	610mm	530mm
Total Length	990mm	1400mm	1230mm	1480mm	1520mm	1550mm	860mm	860mm
Empty Weight	26kg	116kg	116kg	97kg	150kg	165kg	88.5kg	45kg
Max User Weight	114kg	159kg	136kg	182kg	227kg	204kg	136kg	114kg
Battery Capacity	24V/ 11.5A/H	12V/ 75A/H x 2	24V/ 22A/H x 4	24V/ 100A/H x 2	12V/ 75A/H x 2	24V/ 100A/H x 2	12V/ 32A/H x 2	12V/ 14A/H x 2
Turning Diameter	960mm	1855mm	1510mm	1975mm	1650mm	1350mm	500mm	530mm
Ground Clearance	100mm	200mm	120mm	110mm	100mm	130mm	50mm	45mm
Max Safe Slope	6°	10°	10°	10°	8°	10°	6°	6°
Max Driving Range	18km	45.5km	32km	56km	50km	60km	29km	15km
Max Driving Speed	6.4km/h	12km/h	12km/h	12.8km/h	12km/h	15km/h	8km/h	6.4km/h
Motor Capacity	DC 24v 80W Rating 340w max	DC 24v 400W Rating 1750w max	DC 24v 400W Rating 1750w max	DC 24v 650W Rating 2200w max	DC 24v 650W Rating 2200w max	DC 24v 650W Rating 2200w max	DC 24v 160W Rating 520w max	DC 24v 80W Rating 340w max

*Double seat version available

**With lift