



**Goldfern**  
MOBILITY™

# PRODUCT CATALOGUE

# Let us introduce ourselves.

Goldfern is more than a business, it's a commitment to people. Proudly 100% family owned and operated, we are a New Zealand born and bred company based in the Manawatu, providing safe, functional, and affordable solutions that support everyday living.

We understand how important the right equipment and aids are to maintaining independence, comfort, and confidence. Our purpose is simple: to deliver reliable products that integrate seamlessly into daily life. From wheelchairs and mobility scooters to bathroom aids, grab rails, walking aids, and an expanding range of daily living solutions, Goldfern supports people at every stage.

While Goldfern was founded in New Zealand, early visits to our overseas production partners highlighted an opportunity to deliver high-quality products at genuinely accessible prices. Gaining deeper insight into the support and disability sectors revealed rising costs and unmet needs, reinforcing our commitment to do better by listening to end users, healthcare professionals, and resellers.

Every Goldfern product is guided by our uncompromising three-point standard: Safe, Functional, and Affordable. This is more than a principle, it's our promise. While independent safety testing is not a legal requirement, Goldfern is proud to independently test its products against internationally recognised safety standards. It's our way of going above and beyond, giving our customers confidence, peace of mind, and products they can trust.

Backed by strong import and distribution expertise and a growing nationwide network of stockists, Goldfern has built more than a supply chain - we've built trust. We're here to listen, learn, and deliver solutions that genuinely improve everyday life across New Zealand.



Callum Vining  
Director



Michelle Paget  
Key Account Sales Manager

# Table of Contents

|                                        |    |
|----------------------------------------|----|
| <b>Bathroom Aids</b>                   | 04 |
| <b>Daily Living Aids</b>               | 12 |
| <b>Grab Rails</b>                      | 25 |
| <b>Merits Scooter</b>                  | 31 |
| <b>Rollators</b>                       | 43 |
| <b>Walking Aids</b>                    | 54 |
| <b>Wheelchairs</b>                     | 62 |
| <b>Terms of Trade &amp; Warranty's</b> | 71 |

# Bathroom Aids

Designed to enhance independence and safety for individuals with reduced mobility. Our carefully curated selection of bathroom aids ensures a supportive and comfortable environment, promoting confidence and ease of use.

---





## ADJUSTABLE BATH STEP

### FEATURES

- Anodised aluminium frame
- Durable, non-slip plastic platform
- Height adjustable legs
- Non-slip rubber tip/feet

### SPECIFICATIONS

|                   |               |
|-------------------|---------------|
| Adjustable height | 210mm - 260mm |
| Platform size     | 230mm x 410mm |
| Base size         | 285mm x 395mm |
| Max load          | 120kg         |
| Net weight        | 1.15kg        |

SKU GM-YCH-01B04H01



## BATHTUB SEAT

### FEATURES

- Aluminium frame
- Non-slip plastic seat with drainage holes
- Non-slip rubber grips

### SPECIFICATIONS

|                 |               |
|-----------------|---------------|
| Seat size       | 410mm x 230mm |
| Backrest height | 300mm         |
| Overall width   | 740mm         |
| Max load        | 120kg         |
| Net weight      | 1.9kg         |

SKU GM-YCH-01B07H01



## BATH MAT

### FEATURES

- 53x53cm in size
- Perfect for eliminating the risk of slipping
- Attaches to bath surface with suction cups
- Made of PVC material
- Easy to clean

SKU GM-IC-MAT



## FREE STANDING TOILET SURROUND

### FEATURES

- Non-slip feet
- Adjustable width and height
- Padded armrests
- Optional suction feet

### SPECIFICATIONS

|                   |               |
|-------------------|---------------|
| Adjustable height | 675mm - 775mm |
| Adjustable width  | 440mm - 640mm |
| Overall depth     | 490mm         |

### STANDARD

Net weight: 3kg / Max load: 120kg

### HEAVY DUTY

Net weight: 4.2kg / Max load: 200kg



Optional Suction Feet



Adjustable Width

SKU S GM-HMP-8109

SKU HD GM-HMP-8109S

## HEAVY-DUTY BATH BENCH WITH BACK & ARMS

### FEATURES

- Height adjustable design
- Strong cross-brace base
- Durable, non-slip plastic seat and back
- Non-slip feet

**WE'VE MADE SOME CHANGES**  
This product is no longer available



### SPECIFICATIONS

|                             |               |
|-----------------------------|---------------|
| Seat size                   | 305mm x 510mm |
| Width between arms          | 510mm         |
| Net weight                  | 3.55kg        |
| Max load                    | 227kg         |
| Backrest height (from seat) | 335mm         |
| Adjustable height           | 385mm - 515mm |

SKU GM-M-A113

## FOLDING 3-IN-1 COMMODE

### FEATURES

- Durable powder coated steel frame
- Height adjustable legs
- Rubber non-slip feet
- Removable and washable bucket, lid, splash guard and seat
- Removable backrest

### SPECIFICATIONS

|                        |               |
|------------------------|---------------|
| Adjustable seat height | 395mm - 550mm |
| Width between arms     | 455mm         |
| Overall width          | 540mm         |
| Max load               | 113kg         |
| Net weight             | 5.05kg        |

SKU GM-HMP-7210A



## HEAVY-DUTY COMMODE

### FEATURES

- Durable powder coated steel frame
- Height adjustable legs
- Rubber non-slip feet
- Removable and washable bucket with lid and seat
- Simple fold-down arms
- Removable backrest

### SPECIFICATIONS

|                        |               |
|------------------------|---------------|
| Adjustable seat height | 460mm - 570mm |
| Width between arms     | 645mm         |
| Overall width          | 705mm         |
| Max load               | 250kg+        |
| Net weight             | 10.15kg       |

SKU GM-HMP-CCS040






## MOBILE SHOWER COMMUNE

### FEATURES

- Easy-clean powder coated frame
- Moulded seat with removable middle and commode bowl
- Moulded rubber flip-up handles and swing-away leg rests
- Padded backrest
- Locking rear wheels

### SPECIFICATIONS

|                             |                 |
|-----------------------------|-----------------|
| Seat height                 | 530mm           |
| Seat size                   | 440mm x 425mm   |
| Base size                   | 460mm x 840mm   |
| Overall width               | 660mm - 970mm   |
| Armrest height              | 190mm from seat |
| Backrest height             | 415mm from seat |
| Overall height              | 945mm           |
| Net weight (without SP Kit) | 12.85kg         |
| Net weight (with SP Kit)    | 17.6kg          |
| Max load                    | 110kg           |
| Wheels                      | 5" castor       |



Self Propelled Kit  
(Accessory)

SKU GM-AST-MA012-4603

\*Optional Shower Commode Self-Propelled Kit available as requested:

SKU GMP-AST-XLJ-002

## NON-SPILL URINAL BOTTLE

### FEATURES

- Male & female attachments included
- Non-spill design
- Strong plastic
- Autoclavable to 135°
- Easy to clean and hygienic

### SPECIFICATIONS

|          |        |
|----------|--------|
| Capacity | 1000ml |
|----------|--------|

SKU GM-FSE666B



## OVER TOILET FRAME

### FEATURES

- Safe and sturdy support when sitting down and standing up from the toilet
- Lightweight aluminium frame
- Easy clean plastic seat with splash guard
- Height adjustable legs
- Padded armrests

### SPECIFICATIONS

|                    |               |
|--------------------|---------------|
| Adjustable height  | 470mm - 620mm |
| Width between arms | 470mm         |
| Overall width      | 535mm         |
| Net weight         | 3.5kg         |
| Max load           | 120kg         |

SKU GM-HMP7020C



## SHOWER STOOL WITH ROTATING SEAT

### FEATURES

- Round, durable plastic seat
- 360 degree rotating disc on top of seat
- Height adjustable
- Non-slip rubber tip/feet

### SPECIFICATIONS

|                        |               |
|------------------------|---------------|
| Adjustable height      | 340mm x 520mm |
| Outer seat diameter    | 327mm         |
| Rotating disc diameter | 295mm         |
| Max load               | 120kg         |

SKU GM-HMP-5021TF





## SHOWER CHAIR WITH ARMS

### FEATURES

- Epoxy coated mild steel (pre-treated for rust protection) with non-slip rubber tips
- Durable steel arms
- Comfortable non-slip seat and back
- One leg with two-tiered adjustment hole for stability on uneven surfaces
- Drainage holes in seat to prevent water build-up

### SPECIFICATIONS

|                        |               |
|------------------------|---------------|
| Adjustable seat height | 490mm - 640mm |
| Width between arms     | 505mm         |
| Base size              | 535mm x 480mm |
| Max load               | 110kg         |
| Net weight             | 4.5kg         |

SKU GM-HMP792LA



## RAISED TOILET SEAT

### FEATURES

- Easy-clean plastic seat
- Removable cushioned armrests
- Fits most toilets with a secure fastening system
- 3" and 5" option
- Minor assembly required

### SPECIFICATIONS

|          |       |
|----------|-------|
| Max load | 127kg |
|----------|-------|

#### 3" TOILET SEAT

|                    |       |
|--------------------|-------|
| Net weight         | 1.5kg |
| Width between arms | 435mm |

#### 5" TOILET SEAT

|                    |       |
|--------------------|-------|
| Net weight         | 2kg   |
| Width between arms | 450mm |

**NEW SIZES HAVE BEEN ADDED!**  
See product listing for updates

SKU 5" GM-YCH-9002A

SKU 3" GM-YCH-9002



## SHOWER STOOL WITH ARMS

### FEATURES

- Lightweight aluminium
- Durable and strong arms
- Comfortable non-slip seat with drainage
- Height adjustable
- Two-tiered adjustment holes on two legs for stability on uneven surfaces
- Drainage holes in seat to prevent water build-up
- Non-slip rubber tips/feet
- Optional suction feet

### SPECIFICATIONS

|                    |               |
|--------------------|---------------|
| Adjustable height  | 470mm - 610mm |
| Width between arms | 470mm         |
| Seat size          | 285mm x 380mm |
| Max load           | 170kg         |
| Net weight         | 2.3kg         |

SKU GM-HMP792A



**HEAVY-DUTY  
VERSION HAS  
BEEN ADDED!**



## SHOWER STOOL WITHOUT ARMS

### FEATURES

- Lightweight aluminium
- Comfortable non-slip seat with drainage
- Height adjustable
- Two-tiered adjustment holes on two legs for stability on uneven surfaces
- Drainage holes in seat to prevent water build-up
- Non-slip rubber tips/feet
- Optional suction feet

### SPECIFICATIONS

|                   |               |
|-------------------|---------------|
| Adjustable height | 490mm - 610mm |
| Seat size         | 285mm x 380mm |
| Max load          | 170kg         |
| Net weight        | 2kg           |

SKU GM-HMP792

# Daily Living Aids

Our daily living aids are where practicality meets innovation to support individuals with diverse needs in their day-to-day activities. Our curated selection encompasses a wide range of aids designed to enhance independence and convenience for seniors and those with mobility challenges.

---



## BED BACKREST

### FEATURES

- Lightweight frame
- Folds down flat for storage and travel
- Comfortable and easy to clean nylon
- 5 different incline levels

### SPECIFICATIONS

|                  |               |
|------------------|---------------|
| Backrest surface | 570mm x 545mm |
| Base size        | 510mm x 480mm |
| Max load         | 120kg         |
| Net weight       | 3.05kg        |

SKU GM-YCH-8001



Side of Seat

## SWIVEL SEAT

### FEATURES

- Lightweight and easy to transport
- Thick foam cushioning with durable fabric
- Non-slip base to help keep seat in place
- 360° swivel motion

### SPECIFICATIONS

|               |       |
|---------------|-------|
| Seat diameter | 395mm |
| Max load      | 120kg |

SKU GM-YCH-511

### Colours



## JAR OPENER

### FEATURES

- A plastic jar opener for removing pop off lids
- Consists of an integral hook and flat wide handle
- Simple and inexpensive
- Available in red, blue and green

SKU GM-HMP-060



**CURRENTLY UNAVAILABLE**  
New updates coming soon

## BENDABLE COMFORT GRIP CUTLERY

### FEATURES

- Large bendable knife, fork and spoon
- Large bendable handles on all cutlery
- Large bendable heads on fork and spoon
- Made from easy to clean materials
- Comfort-grip handles

SKU F GM-YCH-AWH-88F

SKU K GM-YCH-AWH-88K

SKU S GM-YCH-AWH-88S



**WE'VE MADE SOME CHANGES**  
This product is no longer available

## TABLE TOP JAR OPENER

### FEATURES

- Easy to use with minimal effort
- Holds jars securely in place
- Works on jars with stubborn lids
- Durable plastic construction for everyday use
- Lightweight and easy to handle
- Compact for simple storage and portability

SKU GM-YCH-A5090



## COMB & BRUSH SET

### FEATURES



- Easily cleanable materials
- Soft grip handle
- Soft touch brush bristles

### SPECIFICATIONS

|              |       |
|--------------|-------|
| Comb length  | 375mm |
| Brush length | 350mm |
| Net weight   | 100g  |

SKU GM-IC-8616

## CAR EXIT AID



### FEATURES

- Comfortable contoured handle
- Shaped tip for breaking glass in the event of an emergency
- Hidden blade to assist in cutting the seatbelt in the event of an emergency
- Fits any car door chassis latch
- Stores in the glove box or door pocket

### SPECIFICATIONS

|                |       |
|----------------|-------|
| Max load       | 145kg |
| Overall length | 225mm |

SKU GM-IC-CARAID

## BOTTOM WIPER


### FEATURES

- Easily cleanable materials
- Soft grip handle
- Push button paper release

### SPECIFICATIONS

|                |       |
|----------------|-------|
| Overall length | 395mm |
| Net weight     | 100g  |

SKU GM-IC-8615



## DELUXE REACHER


### FEATURES

- Ergonomic trigger grip with gripped claws, hook, and magnet
- Available in Small or Large size

### SPECIFICATIONS

|                        |       |
|------------------------|-------|
| Overall Length (Large) | 750mm |
| Overall Length (Small) | 600mm |

SKU SMALL GM-HMP-R60  
SKU LARGE GM-HMP-R750





Colours



## DUAL MOTOR ELECTRIC LIFT CHAIR

### FEATURES

- 3M soil and stain inhibitor fabric protection
- Twin motors
- Storage pocket on both sides
- Controller with extra long cord that can be set up for left or right-handed use
- External battery backup
- Stand assist feature

### SPECIFICATIONS

|                            |               |
|----------------------------|---------------|
| Seat depth                 | 530mm         |
| Width between arms**       | 545mm - 700mm |
| Height to seat             | 490mm         |
| Height of back (from seat) | 665mm         |
| Overall width              | 765mm         |
| Overall height             | 1090mm        |
| Base size                  | 615mm x 725mm |
| Max load                   | 150kg         |

\*\*padding and arm bases, not adjustable

SKU CHAIR CHARCOAL

GM-HYE8063

SKU CHAIR DENIM

GM-HYE8063-DEN


SKU CHAIR OATMEAL

GM-HYE8063-OAT

SKU REMOTE


GMP-HYE-8063-L10

### \* LIFT CHAIR PROTECTIVE COVER



Blue

GMP-HYE-COVER-BLUE



Mocha

GMP-HYE-COVER-MOCHA



Grey

GMP-HYE-COVER-GREY



Colours



DAILY LIVING AIDS



## HEAVY-DUTY ELECTRIC LIFT CHAIR

### FEATURES

- New heavy-duty, 3M protected material
- Twin motors
- Storage pocket on both sides
- Controller with extra long cord that can be set up for left or right-handed use
- External battery backup
- User weight of up to 220kg
- Stand assist feature

### SPECIFICATIONS

|                            |               |
|----------------------------|---------------|
| Seat depth                 | 525mm         |
| Width between arms*        | 545mm - 675mm |
| Height to seat             | 460mm         |
| Height of back (from seat) | 645mm         |
| Overall width              | 920mm         |
| Overall height             | 1050mm        |
| Base size                  | 720mm x 920mm |
| Max load                   | 220kg         |

\*padding and arm bases, not adjustable

SKU LIGHT GREY GM-HYE8063HD-LG

SKU STEEL GREY GM-HYE8063HD

Colour



## ZERO GRAVITY ELECTRIC LIFT CHAIR

### FEATURES

- Tilt-in-space function to encourage leg elevation and assist with circulation.
- 3M soil and stain inhibitor fabric protection
- Three motors
- Removable back
- Storage pocket on both sides
- Controller with extra long cord that can be set up for left or right-handed use
- External battery backup
- Stand assist feature
- Available in Charcoal only

### SPECIFICATIONS

|                            |               |
|----------------------------|---------------|
| Seat depth                 | 530mm         |
| Width between arms*        | 545mm - 700mm |
| Height to seat             | 490mm         |
| Height of back (from seat) | 665mm         |
| Overall width              | 765mm         |
| Overall height             | 1090mm        |
| Base size                  | 615mm x 725mm |
| Max load                   | 150kg         |

\*padding and arm bases, not adjustable

SKU GM-HYE8063ZG

## FOLDING PEDAL EXERCISER WITH Pedometer

### FEATURES

- Helps to strengthen arms, hands and legs while increasing blood circulation
- Equipped with a pedometer to track your distance
- Folds down making it easier to slide out of the way when not in use
- Ideal product for use following knee replacement or surgery

### SPECIFICATIONS

|                |                   |
|----------------|-------------------|
| Max load       | 100kg             |
| Overall size   | W 395mm x D 460mm |
| Overall height | 230mm             |
| Net weight     | 2.25kg            |

SKU GM-HMP-S005



## FOOT STOOL WITH HANDLE

### FEATURES

- Strong and durable frame
- Lightweight, easy to manoeuvre
- Non-slip rubber platform surface
- Non-slip tips/ferrules
- Removable handle

### SPECIFICATIONS

|                           |               |
|---------------------------|---------------|
| Platform size             | 360mm x 285mm |
| Base size                 | 345mm x 440mm |
| Platform height           | 230mm         |
| Platform to handle height | 640mm         |
| Floor to handle height    | 870mm         |
| Max load                  | 155kg         |
| Net weight                | 3.65kg        |

SKU GM-YCH-508KD





## GYENNO TWIST SPOON & FORK

### FEATURES

- Self-stabilising mechanics and offsetting abilities for hand tremors up to 70mm
- Support for individuals with hand tremors
- Includes charging cable, portable storage
- Includes fork attachment and spoon attachment

**WE'VE MADE SOME CHANGES**  
This product is no longer available

### SPECIFICATIONS

Net weight 150g

SKU GM-GTS

## INFRARED THERMOMETER

### FEATURES

- Non-contact measurement for hygienic, touch-free use
- Fever alarm - audible alarm if temperature is more than 38°C
- Beep function
- Longevity: 100,000 readings
- Eliminates cross contamination
- Safe for use from newborn onwards
- Three modes: body, room and surface

**WE'VE MADE SOME CHANGES**  
This product is no longer available

### SPECIFICATIONS

|                    |                           |
|--------------------|---------------------------|
| Memory             | 32 readings               |
| Auto power off     | <30 seconds               |
| Measuring distance | 3cm ~ 5cm                 |
| Accuracy           | ± 0.2°C                   |
| Response time      | 1 second                  |
| Power              | 2x AA batteries           |
| Overall size       | H 155mm x W 40mm x D 80mm |
| Net weight         | 80g                       |

SKU GM-GBMD-JXB-178



## LUMBAR SUPPORT CUSHION

### FEATURES

- Contoured shape
- Memory foam
- Perforated breathable fabric
- Adjustable strap

SKU GM-SB-LUMBAR



## SOCK ASSIST

### FEATURES

- Designed to be the no-fuss way to put on socks for those who may otherwise not be able to
- Long fabric handle
- Soft covering makes fitting the socks a no-struggle, comfortable experience

SKU GM-IC-SOCK



## LOTION APPLICATOR

### FEATURES

- Easy way to apply bath gel or body lotion to hard to reach places
- Balls rotate, gently gliding over the skin releasing liquid

### SPECIFICATIONS

Overall length 455mm

SKU GM-IC-LOTION



## OVER WHEELCHAIR TABLE

### FEATURES

- Anodised aluminium frame
- U-shaped base to allow wheelchair access
- Adjustable table top
- Four castor wheels

### SPECIFICATIONS

|                         |                |
|-------------------------|----------------|
| Adjustable height       | 720mm - 1140mm |
| Base size (outside)     | 835mm x 470mm  |
| Base size (inside)      | 735mm x 425mm  |
| Table surface size      | 450mm x 810mm  |
| Table surface thickness | 20mm           |
| Max load                | 40kg           |

SKU GM-HMP-S207

## OVER BED TABLE

### FEATURES

- Split-top configuration: end section remains fixed while larger surface tilts for versatility
- Frame equipped with easy-to-use lift and lower mechanism
- Height adjustable according to preference and bed height
- Ensures ergonomic and comfortable experience for meals, work or reading

### SPECIFICATIONS

|                             |               |
|-----------------------------|---------------|
| Adjustable height           | 710mm - 940mm |
| Overall surface (both tops) | 815mm x 395mm |
| Base size (outside)         | 585mm x 545mm |
| Base size (inside)          | 545mm x 365mm |
| Tilting table surface       | 375mm x 570mm |
| Fixed table surface         | 185mm x 385mm |
| Max load                    | 40kg          |

SKU GM-HMP-S209





## GOLDFERN COMFORT ADULT PULL-UPS

### FEATURES

- Adult pull-ups
- 14 pack
- Unisex
- Wetness indicator
- Odour control
- Latex free and breathable
- Maximum absorbency
- Soft and comfortable


### SPECIFICATIONS

|                                                                                     |      | Waist Size | Est. Absorbency |
|-------------------------------------------------------------------------------------|------|------------|-----------------|
|  | S-M  | 75-100cm   | 1800-1950ml     |
|  | M-L  | 90-120cm   | 1950-2250ml     |
|  | L-XL | 100-150cm  | 2250-2400ml     |

SKU S-M GM-PANTS-S-M

SKU M-L GM-PANTS-M-L

SKU L-XL GM-PANTS-L-XL



## STEPS

### FEATURES

- Rust-free aluminium
- Non-slip mat
- Non-slip rubber feet

### SPECIFICATIONS

|                            |               |
|----------------------------|---------------|
| Foot plate surface         | 450mm x 320mm |
| Height of fixed plate      | 140mm         |
| Height of adjustable plate | 200mm - 250mm |
| Max load                   | 159kg         |
| Net weight                 | 3.65kg        |

|                     |             |
|---------------------|-------------|
| SKU FIXED           | GM-HMP-8334 |
| SKU FIXED W/ HANDLE | GM-HMP-8333 |
| SKU ADJUSTABLE      | GM-HMP-8335 |



## METAL SHOE HORN

### FEATURES

- Extra long shoe horn
- Flexible neck can bend up to 90°
- Designed to be the no-fuss way to put on shoes for those who may otherwise not be able to

SKU GM-YCH-18CS-B



## GRIPMATE

### FEATURES

- 3 piece grip and opening assistance set
- Includes a non-slip mat, multi-size jar opener and multipurpose opener
- Non-slip design for secure handling and better control
- Easy to use and gentle on hands
- Helps reduce effort when opening everyday containers

SKU GM-IC-CR8602-1

# Grab Rails

Goldfern Mobility is dedicated to providing high quality grab rails for those that need a little extra support and stability. We provide a large range of different grab rails for every situation, made of durable materials which are designed to be easy-to-grip.

---





## 45 DEGREE ANGLED GRAB RAIL

### FEATURES

- Stainless Steel or Nylon finish available
- Nylon finish has an antibacterial and anti-uv surface
- Anti-aging surface with texture to add grip
- Strengthened with an aluminium core (Nylon Grab Rails only)
- 45° bend with 2 (310mm) or 3 (600mm) secure mounting bases and 5 fixing positions

### SPECIFICATIONS

Size 600mm x 600mm OR 310mm x 310mm

**STAINLESS STEEL** Circumference: 32mm  
**NYLON** Circumference: 35mm

|                 |                 |
|-----------------|-----------------|
| SKU 600MM SS    | GM-GB45S600X600 |
| SKU 310MM SS    | GM-GB45S310X310 |
| SKU 600MM NYLON | GM-GB45N600X600 |



## FOLD-UP GRAB RAIL

### FEATURES

- Nylon or Stainless Steel finish available
- 90° folding capacity
- Nylon finish is antibacterial, anti-UV and has an anti-aging surface with texture
- Strengthened with an aluminium core (Nylon Grab Rails only)
- Large secure mounting base

### SPECIFICATIONS

#### STAINLESS STEEL

|                |       |
|----------------|-------|
| Height:        | 190mm |
| Width:         | 690mm |
| Circumference: | 32mm  |

#### NYLON

|                |       |
|----------------|-------|
| Height:        | 140mm |
| Width:         | 700mm |
| Circumference: | 35mm  |

|           |               |
|-----------|---------------|
| SKU NYLON | GM-SGBN700    |
| SKU SS    | GM-SGBN700X32 |

## STRAIGHT GRAB RAIL

### FEATURES

- Stainless Steel or Nylon finish available
- Nylon finish is antibacterial, anti-UV and has an anti-aging surface with texture
- Strengthened with aluminium core (Nylon Grab Rails only)
- 2 secure mounting bases with 5 fixing positions
- Available in 8 different sizes



### SPECIFICATIONS

Size 300mm, 450mm, 600mm, 700mm, 800mm, 900mm, 1000mm, 1200mm

|                        |                     |
|------------------------|---------------------|
| <b>STAINLESS STEEL</b> | Circumference: 32mm |
| <b>NYLON</b>           | Circumference: 35mm |

PLEASE VISIT OUR WEBSITE FOR INDIVIDUAL SKU

## L SHAPED GRAB RAIL

### FEATURES

- Nylon or Stainless Steel finish available
- Nylon finish is antibacterial, anti-UV and has an anti-aging surface with texture
- Strengthened with aluminium core (Nylon Grab Rails only)
- 90° bend with 3 secure mounting bases with 5 fixing positions
- Compliant with building code



### SPECIFICATIONS

Size 600mm x 600mm OR 750mm x 750mm

|                        |                     |
|------------------------|---------------------|
| <b>STAINLESS STEEL</b> | Circumference: 32mm |
| <b>NYLON</b>           | Circumference: 35mm |

|                 |                |
|-----------------|----------------|
| SKU 600MM SS    | GM-GBLS600X600 |
| SKU 750MM SS    | GM-GBLS750X750 |
| SKU 600MM NYLON | GM-GBLN600X600 |
| SKU 750MM NYLON | GM-GBLN750X750 |
| SKU 900MM SS    | GM-GBLS900X750 |



Folded Up

## WALL MOUNTED FOLDING SHOWER SEAT 400

### FEATURES

- Quality, antibacterial Nylon surface material
- Lightweight and easy to adjust
- Non-slip seat
- Support leg with non-slip rubber tips
- Fully foldable
- Secured to wall using stainless steel plates

### SPECIFICATIONS

|           |               |
|-----------|---------------|
| Seat size | 400mm x 388mm |
| Max load  | 120kg         |

SKU GM-SSSN400



Folded Up

## WALL MOUNTED FOLDING SHOWER SEAT 845

### FEATURES

- Quality, antibacterial Nylon surface material
- Lightweight and easy to adjust
- Non-slip seat
- Support leg with Nylon tips
- Fully foldable
- Secured to wall using stainless steel plates
- Compliant with building code

### SPECIFICATIONS

|           |               |
|-----------|---------------|
| Seat size | 850mm x 445mm |
| Max load  | 120kg         |

SKU GM-SSSN845



## WALL MOUNTED FOLDING BATH SEAT

### FEATURES

- Quality Nylon surface material
- Nylon finish has an antibacterial and anti-uv surface
- Non-slip seat and tip on support leg
- Lightweight and easy to adjust
- Long lasting stainless steel plates

**WE'VE MADE SOME CHANGES**  
This product is no longer available



Folded Up

### SPECIFICATIONS

|           |               |
|-----------|---------------|
| Seat size | 800mm x 388mm |
| Max load  | 120kg         |

SKU GM-SSSN800



## T-SHAPED GRAB RAIL

### FEATURES

- Stainless Steel available
- Surface texture to add grip and stability
- 3-4 secure mounting bases with 5 fixing positions

### SPECIFICATIONS

|                |       |
|----------------|-------|
| Overall length | 750mm |
| Overall height | 900mm |

**NYLON** Circumference: 35mm

SKU NYLON GM-GBTN900X750



## SHOWER HEAD HOLDER

### FEATURES

- Quality Nylon (strengthened with aluminium core) and Stainless Steel surface materials
- Nylon finish is antibacterial, anti-UV, and anti aging
- Easy to adjust with a wide range of motion
- Fits their corresponding Stainless Steel or Nylon T-Shaped Grab Rail

SKU NYLON    *GM-SHHN*

SKU SS        *GM-SHHSS32*



## FLAT END FOR GRAB RAILS

### FEATURES

- Attaches to all Nylon finish Goldfern Grab Rails
- Used for situations where one end of the grab rail needs to be attached to the adjacent wall

SKU    *GM-FE*

# Merits Scooters

Merits is a brand that stands for faithfulness, care, love and determination. After 30+ years in production, Merits are recognised in various countries around the world as a manufacturer of high quality mobility and rehabilitation products. Merits scooters and power chairs are designed to give users the most comfortable, independent experience possible.

---



Colour



Travel Carry Case

# YOGA PORTABLE TRAVELLING SCOOTER

## FEATURES

- Strong, light-weight alloy aluminium frame
- Automatic safety turning half speed sensor
- Folds down to fit into the included Travel Carry Case
- Height adjustable seat and steering column for maximum user comfort
- Electromagnetic brake
- Available in Red only

## SPECIFICATIONS

|                       |                             |
|-----------------------|-----------------------------|
| Max speed             | 6.4km/h                     |
| Range                 | 18km                        |
| Turn radius           | 960mm                       |
| Ground clearance      | 100mm                       |
| Front wheel           | 6" foam filled tyres        |
| Rear wheel            | 7" foam filled tyres        |
| Motor                 | DC 24V, 80W rating/250W max |
| Battery               | 24V/11.5A/H                 |
| Dimensions            | W 43cm x L 99cm x H 104cm   |
| Overall weight        | 26kg                        |
| Weight in travel case | 31kg                        |
| Max load              | 114kg                       |

SKU    MMS-S542-WR



Colour



MERITSSCOOTERS



Front Section  
Weight - 9.9kg  
(17.4kg with battery)



Rear Section  
Weight - 14.65kg



Seat  
Weight - 12.5kg

## REGAL EZ POWER CHAIR

### FEATURES

- Compact and easy to dismantle for transportation/storage
- Suitable for indoor and limited outdoor use
- High strength steel frame
- Control stick mounting for left and right hand operation
- Thick padded seat
- Intelligent, regenerative electro-magnetic brakes activate automatically when control stick is released
- Adjustable seat, backrest and armrests
- Available in White, Blue, Yellow or Green

### SPECIFICATIONS

|                      |                             |
|----------------------|-----------------------------|
| Max speed            | 6.4km/h                     |
| Range                | 15km                        |
| Turn radius          | 530mm                       |
| Ground clearance     | 45mm                        |
| Wheels               | 200x75 PU tyres             |
| Motor                | DC 24V, 90W rating/300W max |
| Battery              | 12V/14A/H x 2               |
| Dimensions           | W 55cm x L 86cm x H 92cm    |
| Weight incl. battery | 45kg                        |
| Max load             | 114kg                       |

SKU MMS-P321A-WHITE

Colour



## VISION SPORT (LIFT) POWER CHAIR

### FEATURES

- Mid-wheel drive design
- Suitable for indoor and outdoor use
- Full front and rear suspension
- Dual in-line motors
- Adjustable seat, backrest and armrests
- Power seat lift

### SPECIFICATIONS

|                  |                              |
|------------------|------------------------------|
| Max speed        | 8km/h                        |
| Range            | 29km                         |
| Turn radius      | 500mm                        |
| Ground clearance | 50mm                         |
| Wheels           | 10" foam filled tyres        |
| Motor            | DC 24V, 160W rating/520W max |
| Battery          | 12V/32A/H x 2                |
| Dimensions       | W 61cm x L 86cm x H 123cm    |
| Overall weight   | 88.5kg                       |
| Max load         | 136kg                        |

SKU MMS-P326D-WHT



## INDIGO

### FEATURES

- Extendable legroom space of up to 360mm
- Strong lightweight alloy aluminium frame
- Modern LED lighting system
- Digital LCD dashboard display
- Storage compartment under dashboard
- Removable carry basket on front
- Lockable compartment underneath seat
- Electromagnetic brakes
- Wide wheelbase for extra stability
- Adjustable seats, armrests and steering column
- Swivel seat

### SPECIFICATIONS

|                   |                               |
|-------------------|-------------------------------|
| Max speed         | 12km/h                        |
| Range             | 32km                          |
| Turn radius       | 1510mm                        |
| Ground clearance* | Up to 130mm                   |
| Wheels            | 11" tyres                     |
| Motor             | DC 24V, 400W rating/1750W max |
| Battery           | 24V/21A/H x 4                 |
| Dimensions        | W 65cm x L 123cm x H 139cm    |
| Overall weight    | 116kg                         |
| Max load          | 136kg                         |

\*Measured at the highest point of the under-body; additional components on various models may reduce clearance.

SKU GOLD MMS-946A-GOLD

SKU O BLUE MMS-S946A-OB

SKU L BLUE MMS-S946A-LB

SKU D RED MMS-S946A-DR

Colours



## ROADSTER

### FEATURES

- Adjustable seat height and position
- Automatic safety turning half speed sensor
- Removable carry basket with handle
- Generous legroom
- LED lighting system including headlights, tail lights and indicators
- Non-Slip foot mat
- Electromagnetic Brakes
- Height adjustable steering column
- Fully adjustable armrests
- Swivel seat

### SPECIFICATIONS

|                   |                               |
|-------------------|-------------------------------|
| Max speed         | 12km/h                        |
| Range             | 45.5km                        |
| Turn radius       | 1855mm                        |
| Ground clearance* | Up to 200mm                   |
| Wheels            | 12" tyres                     |
| Motor             | DC 24V, 400W rating/1750W max |
| Battery           | 12V/75A/H x 2                 |
| Dimensions        | W 64cm x L 140cm x H 128.5cm  |
| Overall weight    | 115.5kg                       |
| Max load          | 159kg                         |

\*Measured at the highest point of the under-body; additional components on various models may reduce clearance.

SKU GOLD MMS-S745B-G

SKU D BLUE MMS-S745B-DB

SKU D RED MMS-S745B-DR

Colours



MERITS SCOOTERS



## URBAN

### FEATURES

- All round mobility scooter for daily use or longer trips
- Swivel and slide captain's chair for comfort and easy access
- Front and rear suspension for a smooth, stable ride
- Auto headlight and USB charging port
- Supports users up to 182kg
- Up to 56km range with dual batteries
- Adjustable armrests

### SPECIFICATIONS

|                   |                               |
|-------------------|-------------------------------|
| Max speed         | 12.8km/h                      |
| Range             | 56km                          |
| Turn radius       | 1975mm                        |
| Ground clearance* | Up to 150mm                   |
| Wheels            | 14" tyres                     |
| Motor             | DC 24V, 650W rating/2200W max |
| Battery           | 24V/100A/H x 2                |
| Dimensions        | W 66cm x L 148cm x H 124cm    |
| Overall weight    | 158kg                         |
| Max load          | 182kg                         |

\*Measured at the highest point of the under-body; additional components on various models may reduce clearance.

SKU FAIRWAY GREEN MMS-S849A-FG  
SKU LITE BLUE MMS-S849A-LB

SKU TITANIUM MOSS MMS-S849A-TG



Colours



**WE'VE MADE SOME CHANGES**  
This product is no longer available



# SILVERADO EXTREME

## FEATURES

- Sporty streamline design
- Independent front and rear suspension
- Shock absorbers on all four wheels
- Modern LED lighting system
- Digital LCD dashboard display
- Automatic safety turning half speed sensor
- Wide wheelbase for extra stability
- Fully adjustable armrests
- Height adjustable seat and steering column
- Electromagnetic brakes
- Storage compartment under dashboard
- Swivel seat

## SPECIFICATIONS

|                  |                                |
|------------------|--------------------------------|
| Max speed        | 15km/h                         |
| Range            | 60km                           |
| Turn radius      | 1350mm                         |
| Ground clearance | Up to 110mm                    |
| Front wheels     | 14" Pneumatic tyres            |
| Rear wheels      | 16" Pneumatic tyres            |
| Motor            | DC 24V, 650W rating/ 2200W max |
| Battery          | 24V/100A/H x 2                 |
| Dimensions       | W 74cm x L 155cm x H 145cm     |
| Overall weight   | 165kg                          |
| Max load         | 204kg                          |

\*Measured at the highest point of the under-body; additional components on various models may reduce clearance.

SKU D RED    MMS-S940A-DB

SKU SILVER    MMS-S940A-SV



Colours



## CROSS COUNTRY

### FEATURES

- Robust, daring and sporty 4-wheeled outdoor scooter
- Independent front and rear suspension with shock absorbers on all four wheels
- Automatic safety turning half speed sensor
- LED Lighting system
- Electromagnetic brakes
- Storage compartment under dashboard
- Large removable carry basket
- Wide wheelbase for extra stability
- Fully adjustable armrests
- Height adjustable seat and steering column
- Easy-clean removable carpet
- Stylish design
- Swivel seat
- Side-by-side seat "Double Seat" available

### SPECIFICATIONS

|                   |                               |
|-------------------|-------------------------------|
| Max speed         | 12km/h                        |
| Range             | 50km                          |
| Turn radius       | 1650mm                        |
| Ground clearance* | Up to 160mm                   |
| Wheels            | 14" Pneumatic tyres           |
| Motor             | DC 24V, 650W rating/2200W max |
| Battery           | 24V/100A/H x 2                |
| Dimensions        | W 71cm x L 152cm x H 135cm    |
| Overall weight    | 150kg                         |
| Max load          | 227kg                         |

### "DOUBLE SEATER"

Seat base: W 90cm x D 44cm  
 Width between arms: 82.5cm - 91cm  
 Base to top of back: 47cm

\*Measured at the highest point of the under-body; additional components on various models may reduce clearance.

SKU D BLUE MMS-S344A-SB

SKU D RED MMS-S344A-DR

SKU DOUBLE MMS-940-32200385



## ROLLATOR HOLDER

### FEATURES

- Quick and easy fitting via the scooter's rear attachment slot
- Holds all current Goldfern rollator models

SKU GM-MMS-RH



## SCOOTER BACKPACK WITH FRAME

### FEATURES

- Waterproof material with reflective rear safety strip
- Velcro fastening for quick and easy access
- Highly durable, adjustable frame kit
- Compatible with all scooter models
- Secure cane holder (cane not included)

SKU MMS-34200297



## REAR FLAG

### FEATURES

- Orange, bright-coloured flag to add visibility to mobility scooters
- Fits easily to all scooters via the rear attachment slot

CONTACT US DIRECTLY TO ORDER



## THE CRUISER

### FEATURES

- Can fit in the back of an SUV or station-wagon to transport with ease
- Optional umbrella
- Optional rear seat with storage
- Front golf bag holder
- Electromagnetic braking
- Charges in only 4 hours
- Easy-to-learn comfortable and stable steering mechanics
- 80% more gentle on the green compared to the traditional golf cart

### SPECIFICATIONS

|                     |                 |
|---------------------|-----------------|
| Dimensions          | 1280x720x1300mm |
| Speed options       | 15km/h & 20km/h |
| Max load            | 150kg           |
| Max range           | 30+km           |
| Turn radius         | 120cm           |
| Ground clearance    | 110mm           |
| Suspension          | Centre shock    |
| Tyres               | 10 inch         |
| Working temperature | -10°C to 40°C   |

CONTACT US DIRECTLY TO ORDER



Premium Cart



Standard Cart



Cart with Tray

## EURO GOLF CARTS

### MODELS

- Standard Cart
- Premium Cart
- Standard Cart with Tray

### PREMIUM CART FEATURES

- Ball Washer
- Chilly Bin
- Luxury Seats
- Headlights
- Rain Cover
- Black Roof

### SPECIFICATIONS

|                       |                   |
|-----------------------|-------------------|
| Dimensions            | 2400x1245x1900mm  |
| Motor                 | AC Motor 48v 5kw  |
| Battery               | Lithium Battery   |
| Charger               | Tonsload Onboard  |
| Net weight            | 430kg             |
| Climbing capacity     | 30%               |
| Ground clearance      | 114mm             |
| Payload               | 360kg             |
| Max speed             | 25km/h            |
| Cart tray material    | Aluminium         |
| Cart tray dimensions* | L 1.73m x W 1.15m |
| Cart tray depth       | 255mm             |

CONTACT US DIRECTLY TO ORDER

\*Tray size can be customised on request

# Goldfern Rollators

Designed to empower individuals with stability, support and style. Our range showcases a variety of rollators engineered to meet diverse needs and preferences. Whether you're seeking lightweight models for easy maneuverability or robust designs for outdoor adventures, our range offers unparalleled comfort and convenience. With features such as adjustable height, ergonomic handles and secure braking systems, our rollators prioritise safety without compromising on style.

---



Please note: product weight may vary due to manufacturing differences and optional accessories.

Colour



Folded Down

## 3 WHEELED ALUMINIUM ROLLATOR

### FEATURES

- Aluminium tube frame
- Ergonomic, height adjustable handles with locking handbrakes
- Swivel front wheel
- Vinyl storage bag with zip closure
- Easily manoeuvres in tight spaces
- Folds up for transport, storage and standby
- Suitable for indoor use

### SPECIFICATIONS

|                          |               |
|--------------------------|---------------|
| Adjustable handle height | 820mm - 945mm |
| Wheels                   | 8"            |
| Base width               | 620mm         |
| Base depth               | 580mm         |
| Net weight               | 4.65kg        |
| Max load                 | 114kg         |

SKU GM-M-W340-2




Colours



Curb Climber



Fasten frame shut



Clip frame open

ROLLATORS

## ECO-LITE ROLLATOR

### FEATURES

- Aluminium frame
- Curved folding backrest
- Nylon under-seat bag
- Walking stick holder
- Optional food tray
- Curb climber
- Locking clip to secure when open
- Strap to secure shut when folded
- Soft compound ergonomic handles
- Long Handle Kit available as accessory
- ISO 11199-2:2005 tested and SGS approved
- Suitable for indoor and outdoor use
- Long handle kit accessory available

### SPECIFICATIONS

|                          |               |
|--------------------------|---------------|
| Adjustable handle height | 845mm - 935mm |
| Width between handles    | 430mm         |
| Overall width            | 580mm         |
| Overall depth            | 695mm         |
| Seat height              | 560mm         |
| Wheels                   | 8"            |
| Net weight               | 6.3kg         |
| Max load                 | 136kg         |

SKU RED GM-TMA-8050ELR  
SKU BLUE GM-TMA-8050ELB

SKU GREEN GM-TMA-8050ELG

SKU PURPLE GM-TMA-8050ELP



Folded Down



**WE'VE UPDATED OUR BAGS!**  
New bag shown on blue version

Colours



## X-FOLD™ ROLLATOR

### FEATURES

- Suitable for indoor and outdoor use
- Removable shopping bag with easy to open zipper
- Lightweight aluminium frame
- Secure opening/closing mechanisms
- Removable wheels for travelling
- Walking stick holder
- Curb climber

### SPECIFICATIONS

|                            |                |
|----------------------------|----------------|
| Adjustable handle height   | 895mm - 1020mm |
| Width between handles      | 460mm          |
| Overall width              | 610mm          |
| Overall depth              | 735mm          |
| Seat size                  | 455mm x 200mm  |
| Seat height                | 560mm          |
| Wheels                     | 8"             |
| Net weight excl. Carry Bag | 5.85kg         |
| Max load                   | 136kg          |

SKU GREEN GM-TMA-KR56-GREEN  
 SKU GREY GM-TMA-KR56-GREY  
 SKU ICE BLUE GM-TMA-KR56-ICEBLUE

SKU PURPLE GM-TMA-KR56-PURPLE  
 SKU PINK GM-TMA-KR56-PINK



Colours



ROLLATORS

Folded Down

## X-FOLD™ MINI ROLLATOR

### FEATURES

- Suitable for indoor and outdoor use
- Removable shopping bag with shoulder strap for easy carrying
- Lightweight aluminium frame
- Secure opening/closing mechanisms
- Removable wheels for travelling
- Walking stick holder
- Curb climber

### SPECIFICATIONS

|                            |               |
|----------------------------|---------------|
| Adjustable handle height   | 805mm - 880mm |
| Width between handles      | 385mm         |
| Overall width              | 535mm         |
| Overall depth              | 725mm         |
| Seat size                  | 370mm x 200mm |
| Seat height                | 515mm         |
| Wheels                     | 8"            |
| Net weight excl. Carry Bag | 5.65kg        |
| Max load                   | 136kg         |

SKU GREEN GM-TMA-KR51N-GREEN  
SKU GREY GM-TMA-KR51N-GREY

SKU PURPLE GM-TMA-KR51N-PURPLE

Colour



## X-FOLD™ PLUS ROLLATOR

### FEATURES

- Suitable for indoor and outdoor use
- Removable shopping bag with shoulder strap for easy carrying
- High-strength aluminium frame design
- Secure opening/closing mechanisms
- Removable wheels for travelling
- Walking stick holder
- Curb climber

### SPECIFICATIONS

|                            |               |
|----------------------------|---------------|
| Adjustable handle height   | 810 - 935mm   |
| Width between handles      | 525mm         |
| Overall width              | 655mm         |
| Overall depth              | 725mm         |
| Seat size                  | 500mm x 200mm |
| Seat height                | 560mm         |
| Wheels                     | 8"            |
| Net weight excl. Carry Bag | 6.3kg         |
| Max load                   | 180kg         |

SKU GM-TMA-KR56W



Colour



ROLLATORS

## HEAVY-DUTY PREMIUM ROLLATOR

### FEATURES

- Fully foldable powder coated steel frame
- Wide rear frame to maximise stride space
- Lockable rear mounted hand brakes
- Suited to indoor and outdoor use
- Adjustable handle height
- Comfortable foam padded backrest
- Vinyl under seat bag
- Optional food tray
- Long handle kit accessory available

### SPECIFICATIONS

|                          |               |
|--------------------------|---------------|
| Adjustable handle height | 855mm - 942mm |
| Width between handles    | 490mm         |
| Overall width            | 635mm         |
| Overall depth            | 690mm         |
| Seat height              | 580mm         |
| Wheels                   | 8"            |
| Net weight               | 10.8kg        |
| Max load                 | 180kg         |

SKU GM-TMA8350RIT-HD



Colour



## PREMIUM 8" ROLLATOR

### FEATURES

- Foldable lightweight aluminium frame
- Handbrake locking mechanism
- Adjustable handle height
- Padded back bar
- Vinyl under seat bag
- Long Handle Kit available as accessory
- ISO 11199-2 tested and approved
- Optional food tray
- Available with 6" or 8" wheels (6" wheels for indoor use only)
- Long handle kit accessory available

### SPECIFICATIONS

|                          |               |
|--------------------------|---------------|
| Adjustable handle height | 855mm - 940mm |
| Width between handles    | 440mm         |
| Overall width            | 560mm         |
| Overall depth            | 685mm         |
| Seat size                | 340mm x 320mm |
| Seat height              | 570mm         |
| Wheels                   | 6"/8"         |
| Net weight               | 7.65kg        |
| Max load                 | 127kg         |

SKU 6" GM-TMA8350RITBLK

SKU 8" GM-TMA8350RIT185



Colour



ROLLATORS

## PREMIUM NARROW 8" ROLLATOR

### FEATURES

- Lightweight foldable aluminium frame
- Folding frame design for storage/transport
- 8" Polyurethane wheels suitable for both indoor and outdoor use
- Height-adjustable ergonomic handles
- Handbrakes with rear wheel parking lock
- Vinyl under seat storage bag
- Optional food tray and cane holder
- ISO 11199-2:2021 tested and approved

### SPECIFICATIONS

|                          |               |
|--------------------------|---------------|
| Adjustable handle height | 800mm - 890mm |
| Width between handles    | 380mm         |
| Overall width            | 550mm         |
| Overall depth            | 670mm         |
| Seat size                | 340mm x 305mm |
| Seat height              | 510mm         |
| Wheels                   | 8"            |
| Net weight               | 7.3kg         |
| Max load                 | 127kg         |

SKU GM-TMA8350RITN185



In use



Folded down

## UNIVERSAL FOOD TRAY

### FEATURES

- Fits the Premium Rollator range (GM-TMA8350RIT/RIT185/RIT185HD/RITN185) and the ECO-Lite Rollator (GM-TMA-8050)

### SPECIFICATIONS

|                |               |
|----------------|---------------|
| Tray size      | 405mm x 280mm |
| Tray thickness | 30mm          |

AVAILABLE ON REQUEST



## LONG HANDLE KIT

### FEATURES

- Designed to fit the X-Fold walker range
- Longer handles reduce bending and improve posture
- Provides a more natural and comfortable walking position
- Ideal for taller users
- Quick and easy to install without replacing the full walker
- Height adjustable for a personalised fit

### SPECIFICATIONS

|                          |                |
|--------------------------|----------------|
| Adjustable handle height | 890mm - 1010mm |
|--------------------------|----------------|

### COMPATIBLE WITH:

X-Fold Rollator  
X-Fold PLUS Rollator

**\*Kit only, rollators sold separately**

SKU GMP-TMALHKKR56



## ROLLATOR CANE HOLDER

### FEATURES

- Attaches to all standard A-Frame rollators

SKU GMP-TMA-18-008



## GUTTER ARM KIT

### FEATURES

- Designed to fit the X-Fold walker range
- Allows better weight distribution to be supported through the forearms, reducing strain on hands and wrists
- Ideal for users with arthritis, weak grip, or hip/ pelvic recovery
- Promotes a more upright posture and better balance
- Includes arm pads/brake/post, brake cable and rear wheels
- Easy to fit with some assembly required

### SPECIFICATIONS

Adjustable handle height 1080mm - 1180mm

*\*When fitted to a standard X-Fold Rollator*

Handle reach 265mm - 365mm

*\*Taken from end of arm pad to brake handle*

### COMPATIBLE WITH:

X-Fold Rollator

X-Fold PLUS Rollator

**\*Kit only, rollators sold separately**

SKU GMP-TMA-23S-035



# Walking Aids

Experience newfound confidence and support in every step with our selection of rollators, canes, and more. Our thoughtfully curated range provides stability, comfort, and style, allowing individuals to maintain an active and independent lifestyle. Whether strolling through the neighbourhood or tackling daily tasks, our walking aids offer the perfect blend of functionality and design to enhance mobility and promote a life in motion.

---





## CURVED CROOK HANDLE WALKING STICK

### FEATURES

- Lightweight aluminium
- Crook shaped handle with foam grip
- 10 adjustable height levels with lock
- Non-slip rubber tip

### SPECIFICATIONS

|                   |                            |
|-------------------|----------------------------|
| Adjustable height | 812mm - 914mm              |
| Tip diameter      | Inner - 19mm / Base - 40mm |
| Max load          | 113kg                      |

SKU BLACK GM-HMP-20720BLK

SKU BRONZE GM-HMP-20720BRZ



## CURVED SWAN HANDLE WALKING STICK

### FEATURES

- Lightweight aluminium
- Curved swan shaped handle with foam grip
- 10 adjustable height levels with lock
- Non-slip rubber tip
- Wrist strap on handle end

### SPECIFICATIONS

|                   |                            |
|-------------------|----------------------------|
| Adjustable height | 787mm - 1016mm             |
| Tip diameter      | Inner - 19mm / Base - 40mm |
| Max load          | 113kg                      |

SKU LEOPARD GM-HMP-20042LEO

SKU BLACK GM-HMP-20042BLK

SKU BRONZE GM-HMP-20042BRZ



KIWIANA

BLACK

HUE TETOKO

FRESHIA

BRONZE

RED BONNET

CHEETAH

## FOLDING WALKING STICK

### FEATURES

- Lightweight aluminium
- Comfortable wooden T shaped handle
- 10 adjustable height levels with lock
- Non-slip rubber tip
- Wrist strap below handle

### SPECIFICATIONS

|                   |                            |
|-------------------|----------------------------|
| Adjustable height | 812mm - 914mm              |
| Tip diameter      | Inner - 19mm / Base - 40mm |
| Max load          | 113kg                      |

|                |                   |
|----------------|-------------------|
| SKU KIWIANA    | GM-HMP-20061KIWI  |
| SKU BLACK      | GM-HMP-20061BLK   |
| SKU HUE TETOKO | GM-HMP-20061HUE   |
| SKU CHEETAH    | GM-HMP-20061CHEET |

|                |                  |
|----------------|------------------|
| SKU FRESHIA    | GM-HMP-20061-193 |
| SKU BRONZE     | GM-HMP-20061BRZ  |
| SKU RED BONNET | GM-HMP-20061RH   |



## SELF STANDING BASE

### FEATURES

- Allows stick to stand upright when not in use
- Designed to fit 19mm canes
- Made of flexible black rubber

SKU GM-IC-SELFSTANDBASE

## T HANDLE WALKING STICK

### FEATURES

- Lightweight aluminium
- Comfortable wooden T shaped handle
- 10 adjustable height levels with lock
- Non-slip rubber tip
- Wrist strap below handle

### SPECIFICATIONS

|                   |                            |
|-------------------|----------------------------|
| Adjustable height | 736mm - 965mm              |
| Tip diameter      | Inner - 19mm / Base - 40mm |
| Max load          | 113kg                      |



|            |                   |
|------------|-------------------|
| SKU BLOOM  | GM-HMP-20052-P081 |
| SKU PURPLE | GM-HMP-20052-151  |
| SKU BRONZE | GM-HMP-20052BRZ   |
| SKU BLACK  | GM-HMP-20052BLK   |
| SKU ROSE   | GM-HMP-20052-ROSE |

## QUAD STICK SMALL & LARGE

### FEATURES

- Lightweight aluminium with swan neck
- Comfortable foam handle
- 10 options for adjustable height
- Non-slip rubber tip/ferrule
- Available with small or large base size

### SPECIFICATIONS

|                   |                            |
|-------------------|----------------------------|
| Adjustable height | 736mm - 965mm              |
| Tip diameter      | Inner - 13mm / Base - 25mm |
| Net weight        | 800g                       |
| Max load          | 100kg                      |

|              |                          |
|--------------|--------------------------|
| <b>SMALL</b> | Base Size: 152mm x 203mm |
| <b>LARGE</b> | Base Size: 190mm x 304mm |



|           |              |
|-----------|--------------|
| SKU SMALL | GM-HMP-40022 |
| SKU LARGE | GM-HMP-40021 |



## HEAVY-DUTY WALKING STICK

### FEATURES

- Soft cushioned handle with wrist strap
- Adjustable height with dual locking
- Non-slip rubber tip/ferrule

### SPECIFICATIONS

|                   |                          |
|-------------------|--------------------------|
| Adjustable height | 740mm - 1020mm           |
| Tip diameter      | Inner - 22 / Base - 50mm |
| Max load          | 250kg                    |

SKU GM-HMP-20045BLK



Folded down

## TRIPOD WALKING STICK WITH SEAT

### FEATURES

- Adjustable height
- Lightweight aluminium frame
- Folding design opens up to create a sturdy seat
- Comfortable plastic moulded seat design

### SPECIFICATIONS

|                   |               |
|-------------------|---------------|
| Adjustable height | 860mm - 990mm |
| Seat height       | 505mm - 605mm |
| Net weight        | 1.3kg         |
| Max load          | 100kg         |

SKU GM-HMP-21001



## TABLE TOP CANE HOLDER

### FEATURES

- Clips on to most standard walking sticks
- Folds easily
- Perches perfectly on table or solid surface edges

SKU GM-IC-CH002



## ADULT STANDARD CRUTCH

### FEATURES

- Height-adjustable handgrip and round cuff
- ISO 11334-1:2007 tested and approved
- Sold individually

### SPECIFICATIONS

|                         |               |
|-------------------------|---------------|
| Base to handgrip height | 710mm - 940mm |
| Handgrip to cuff height | 255mm - 325mm |

SKU GM-HMP-3005AN



## JUNIOR GRIP CRUTCH

### FEATURES

- Height-adjustable handgrip and oval cuff
- Cuff design tailored for youth users
- Sold individually

### SPECIFICATIONS

|                         |               |
|-------------------------|---------------|
| Base to handgrip height | 550mm - 780mm |
| Handgrip to cuff height | 250mm - 320mm |

SKU GM-HMP-3005YN



## ADULT ANATOMICAL CRUTCHES

### FEATURES

- Height-adjustable handgrip and cuff
- Ergonomic moulded plastic handgrip
- ISO 11334-1:2007 tested and approved
- Sold as a pair

### SPECIFICATIONS

|                         |                |
|-------------------------|----------------|
| Base to handgrip height | 795mm - 1010mm |
| Handgrip to cuff height | 280mm - 355mm  |

SKU GM-HMP-3008A

Colour




Folded Down

## STANDARD KNEE WALKER

### FEATURES

- Adjustable handle height
- Suitable for indoor and outdoor use
- Dual-locking hand brakes
- Easy to fold down handle (fits below seat)
- Front carry basket
- Easy to adjust knee pad height
- Heavy Duty version available

### SPECIFICATIONS

|                             |                |
|-----------------------------|----------------|
| Adjustable handle height    | 825mm - 1060mm |
| Overall width               | 515mm          |
| Wheel to wheel width        | 400mm          |
| Overall length excl. basket | 800mm          |
| Knee pad height             | 370mm - 540mm  |
| Knee pad size               | 185mm x 370mm  |
| Wheels                      | 8"             |
| Net weight incl. basket     | 11.3kg         |
| Max load                    | 159kg          |

SKU GM-UN-3023

Colour



Folded Down

## HEAVY-DUTY KNEE WALKER

### FEATURES

- Comfortable anatomical hand grips
- Adjustable handle height
- Suitable for indoor and outdoor use
- Dual-locking hand brakes
- Easy to fold down handle (fits below seat)
- Front carry basket
- Easy to adjust knee pad height
- Reflective rear sticker
- Standard version available

### SPECIFICATIONS

|                          |                |
|--------------------------|----------------|
| Adjustable handle height | 825mm - 1060mm |
| Overall width            | 400mm          |
| Overall length           | 800mm          |
| Knee pad height          | 400mm - 560mm  |
| Knee pad size            | 185mm x 370mm  |
| Wheels                   | 8"             |
| Net weight incl. basket  | 11.6kg         |
| Max load                 | 180kg          |

SKU GM-UN-3024

# Goldfern Wheelchairs

Experience freedom without limits with our collection of wheelchairs, meticulously crafted to redefine mobility and accessibility. From lightweight transport chairs ideal for travel, to self-propelled heavy-duty wheelchairs tailored to unique needs. Our catalogue offers a comprehensive range of options designed to enhance comfort, maneuverability and independence.

---





## ECONOMY TRANSIT WHEELCHAIR

### FEATURES

- Steel folding frame
- Black nylon upholstery
- Fixed padded armrests
- Removable swing-away footrests

### SPECIFICATIONS

|                         |                             |
|-------------------------|-----------------------------|
| Dimensions*             | W 605mm x H 910mm x D 730mm |
| Front wheels            | 6" solid rubber             |
| Rear wheels             | 8" solid rubber             |
| Back height (from seat) | 410mm                       |
| Seat size               | 18" x 16"                   |
| Floor to seat height    | 440mm                       |
| Net weight*             | 10.2kg                      |
| Max load                | 136kg                       |

\*Excluding footrests

SKU GM-AST-MS019-003-1

Please note that dimensions may vary by up to 5mm



Folded Up

## ECONOMY SELF-PROPELLED WHEELCHAIR

### FEATURES

- Steel folding frame
- Black nylon upholstery
- Fixed padded armrests
- Fixed and fold-up footrests

### SPECIFICATIONS

|                         |                             |
|-------------------------|-----------------------------|
| Dimensions              | W 640mm x H 915mm x D 790mm |
| Front wheels            | 6" solid rubber             |
| Rear wheels             | 24" solid rubber            |
| Back height (from seat) | 430mm                       |
| Seat size               | 18" x 16"                   |
| Floor to seat height    | 440mm                       |
| Net weight              | 15.3kg                      |
| Max load                | 136kg                       |

SKU GM-AST-MS001-SP

Please note that dimensions may vary by up to 5mm



## LIGHTWEIGHT TRANSIT WHEELCHAIR

### FEATURES

- Removable swing-away footrests
- Height adjustable foot plates and arm pads
- Adjustable anti-tippers and lap belt
- Folding frame
- Quick release wheels
- Tested and approved to AS/NZS 3695.1:2011

### SPECIFICATIONS

|                         |                             |
|-------------------------|-----------------------------|
| Dimensions*             | W 645mm x H 870mm x D 735mm |
| Front wheels            | 6" solid rubber             |
| Rear wheels             | 12" solid rubber            |
| Back height (from seat) | 380mm                       |
| Seat size               | 18" x 16"                   |
| Width between armrests  | 460mm                       |
| Floor to seat height    | 470mm                       |
| Net weight*             | 9.45kg                      |
| Max load                | 100kg                       |

\*Excluding footrest

SKU GM-AST-LIGHT-AT

Please note that dimensions may vary by up to 5mm



## LIGHTWEIGHT SELF-PROPELLED WHEELCHAIR

### FEATURES

- Aluminium folding frame
- Folding nylon padded backrest and seat
- Attendant brakes and user wheel brakes
- Removable swing-away footrests
- Lap belt and rear anti-tippers
- Quick release wheels
- Tested and approved to AS/NZS 3695.1

### SPECIFICATIONS

|                         |                             |
|-------------------------|-----------------------------|
| Dimensions*             | W 680mm x H 870mm x D 735mm |
| Front wheels            | 6" solid rubber             |
| Rear wheels             | 20" quick release           |
| Back height (from seat) | 380mm                       |
| Width between arms      | 460mm                       |
| Seat size               | 18" x 16"                   |
| Floor to seat height    | 470mm                       |
| Net weight*             | 10.1kg                      |
| Max load                | 100kg                       |

\*Excluding footrest

SKU GM-AST-LIGHT-ASP

Please note that dimensions may vary by up to 5mm



## PREMIUM TRANSIT WHEELCHAIR

### FEATURES

- Heavy duty nylon upholstery
- Adjustable anti-tip bars to prevent wheelchair tipping
- Lap belt
- Main frame folds inwards
- Detachable desk armrests
- Height and length adjustable arm pads
- Removable swing-away footrests
- Adjustable and removable calf strap
- Fixed backrest handles
- Solid rubber wheels (8" front, 12" back)
- Height adjustable seat
- Optional elevating leg rests and stump support
- Available in 2 seat sizes

### SPECIFICATIONS

|                        |         |
|------------------------|---------|
| Net weigh*             | 11.80kg |
| Max load               | 125kg   |
| <b>16"</b>             |         |
| Depth*                 | 775mm   |
| Width                  | 580mm   |
| Floor to handle height | 915mm   |
| Floor to seat height   | 515mm   |
| Width between arms     | 415mm   |
| <b>18"</b>             |         |
| Depth*                 | 780mm   |
| Width                  | 605mm   |
| Floor to handle height | 905mm   |
| Floor to seat height   | 515mm   |
| Width between arms     | 465mm   |

\*Excluding footrest

SKU 16" GME-AST-EUT1616

SKU 18" GME-AST-EUT1818

Please note that dimensions may vary by up to 5mm



## PREMIUM SELF-PROPELLED WHEELCHAIR

### FEATURES

- Heavy duty nylon upholstery
- Adjustable anti-tip bars to prevent wheelchair tipping
- Lap belt
- Detachable desk armrests
- Height and length adjustable arm pads
- Removable swing-away footrests
- Adjustable and removable calf strap
- Fixed backrest handles
- Solid rubber wheels (8" front, 24" back)
- Height adjustable seat
- Optional elevating leg rests and stump support
- Quick release wheels
- Tested and approved to AS/NZS 3695.1:2011
- Available in 2 seat sizes

### SPECIFICATIONS

|                        |        |
|------------------------|--------|
| Net weight*            | 13.5kg |
| Max load               | 125kg  |
| <b>16"</b>             |        |
| Depth*                 | 800mm  |
| Width                  | 610mm  |
| Floor to handle height | 930mm  |
| Floor to seat height   | 515mm  |
| Width between arms     | 415mm  |
| <b>18"</b>             |        |
| Depth*                 | 800mm  |
| Width                  | 660mm  |
| Floor to handle height | 930mm  |
| Floor to seat height   | 515mm  |
| Width between arms     | 460mm  |

\*Excluding footrest

SKU 16" GME-AST-EUSP1616

SKU 18" GME-AST-EUSP1818

Please note that dimensions may vary by up to 5mm



## HEAVY-DUTY WHEELCHAIR

### FEATURES

- Double cross-braces for extra strength, rigidity and weight capacity
- Dual axles provide semi-height convertibility, allowing seat re-positioning
- Detachable desk armrests or full armrests
- Removable swing-away footrests made of aluminium
- Nylon seat and back upholstery

### SPECIFICATIONS

|                         |                                |
|-------------------------|--------------------------------|
| Dimensions              | W 560mm x L 1070mm x D 835mm   |
| Adjustable height       | 840mm - 890mm                  |
| Width between arms      | 565mm                          |
| Front wheels            | 8" x 1-1/4" solid tyre casters |
| Rear wheels             | 24" x 1" solid tyres mag wheel |
| Back height (from seat) | 390mm                          |
| Seat size               | 565mm x 460mm                  |
| Floor to seat height    | 470mm - 520mm                  |
| Net weight              | 24.6kg                         |
| Max load                | 227kg                          |

SKU GM-M-N473-22"

Please note that dimensions may vary by up to 5mm



## PAEDIATRIC WHEELCHAIR

### FEATURES

- Removable swing-away footrests
- Height adjustable footrests and arm pads
- Adjustable anti-tippers and lap belt
- Height adjustable attendant handles
- Folding frame
- Quick release wheels
- Tested and approved to AS/NZS 3695.1:2011
- Available in 2 seat sizes

### SPECIFICATIONS

|                           |                             |
|---------------------------|-----------------------------|
| Dimensions                | W 485mm x H 800mm x D 710mm |
| Width between arms        | 315mm                       |
| Adjustable handles        | 800mm - 940mm               |
| Front wheels              | 6" solid rubber             |
| Rear wheels               | 20" quick-release           |
| Back height (from seat)   | 360mm                       |
| Seat size                 | 12" x 14" OR 14" x 14"      |
| Floor to seat height      | 395mm                       |
| Net weight excl. footrest | 12.5kg                      |
| Max load                  | 100kg                       |

SKU 12" GM-AST-JNRSP-1212  
SKU 14" GM-AST-JNRSP-1414

Please note that dimensions may vary by up to 5mm

# TERMS OF TRADE, WARRANTY'S & WORKMANSHIP GUIDE

---

## DELIVERY

Goldfern Mobility work closely with two major freight carriers. All of the products featured here are available via our Palmerston North based warehouse and can be delivered nationwide. Contact one of our team for shipping costs. Stockists will enjoy FREE FREIGHT on all orders over the minimum spend amount\*.

For all orders under the minimum spend amount, a flat rate freight charge will be applied.\*

*\* These freight charges are subject to change, please see our website for up to date freight charges and minimum spend amounts.*



### INTERNATIONAL CUSTOMERS

Currently we are only shipping to New Zealand addresses.



### RETURNS

Please choose carefully as Goldfern Mobility cannot accept goods for return if you have changed your mind. For any goods found to be faulty we will fulfill our obligations under the Consumers Guarantees Act 1993.



### SAFETY STANDARDS

Our goal is to be affordable and functional, but we consider the safety of those using our product to be paramount.

## ISO & AS/NZS APPROVED PRODUCTS

Should you see an item with one of these stamps, you can be sure it has been tested and approved by an SGS testing lab to strict, internationally recognised production and safety standards.

Curious about these standards, testing procedures and safety of our products? We would love for you to get in touch to discuss our products with you further.



---

## PARTS AVAILABLE

We stock a wide selection of replacement parts and accessories to support our products long after purchase. This helps ensure your equipment can be maintained, repaired, and kept in good working order, providing confidence in long-term use and reliability.

# WARRANTY POLICY

Goldfern Mobility provides a warranty on most of our products from date of purchase. Please refer to each product specification for the warranty period applicable. This warranty is to provide a guarantee that any product purchased from Goldfern Mobility will be free from material and workmanship defects.

If a problem does arise with any Goldfern Mobility product within the warranty period, please complete the warranty/repair form on our website including images. The form can be downloaded from our website [www.goldfernmobility.co.nz](http://www.goldfernmobility.co.nz).

Alternatively, you can email us at [hello@goldfernmobility.co.nz](mailto:hello@goldfernmobility.co.nz), or phone us on **0800 100 971** to receive a return authorisation.

Depending on the circumstance and product, Goldfern Mobility may request for the product to be returned to our Head Office for assessment or we may arrange for a sub-contractor to inspect on site.

Should a product be found to be defective, Goldfern Mobility reserves the right to either repair or replace the item. We endeavour to resolve all warranty issues as quickly as possible.

Warranties will not apply to products that have been damaged due to normal wear and tear, when products have not been used as intended, when unapproved modifications have been made to products, or been used for hire.

Goldfern Mobility agree to fulfil our obligations under the Consumer Guarantees Act 1993.

Products used in any form of hire pool will be subject to relative warranty terms and conditions. Please discuss this with Goldfern Mobility.

**For Product Warranty's please contact us or refer to our website: [www.goldfernmobility.co.nz](http://www.goldfernmobility.co.nz)**

**Items in this catalogue are subject to change - please refer to our website, or call us prior to confirming your order.**

## SCOOTER WARRANTY

Goldfern Mobility, on behalf of Merits Health Products, warrants these Mobility Scooters and Power Chairs against defects in materials and workmanship for a designated period from the date of purchase. Warranties are subject to certain limitations and exclusions as defined below and supported throughout New Zealand by Goldfern Mobility Ltd.

Goldfern Mobility will, at their sole discretion, repair or replace parts found upon examination by an authorised representative of Goldfern Mobility, to be defective in material and/or workmanship.

### TWO YEAR LIMITED WARRANTY

Structural frame components, including:

- Main frame
- Seat post
- Fork
- Tiller frame

### ONE YEAR LIMITED WARRANTY

Drive components, including:

- Motors
- Transaxle
- Brake
- Electronics

Battery and Battery Charger components supplied with your Merits Mobility Scooter or Power Chair are covered by a separate warranty offered by the battery and battery charger manufacturer. This is for a period outlined below from the date of purchase by the reseller. In addition, Goldfern Mobility will provide an extended pro-rata warranty for a further 6 months, providing cover in the event of failure as follows:

#### 0-12 MONTHS

100% replacement of parts cost

#### 12-15 MONTHS

67% replacement of parts cost

#### 15-18 MONTHS

33% replacement of parts cost

# GRAB RAIL WARRANTY

This warranty covers items within the Grab Rail range, including nylon or stainless-steel rails and nylon bath/shower seats.

The product or product component is covered by a 12 month warranty from date of purchase by the reseller for faulty workmanship or materials subject to correct installation by a qualified plumber, builder or relevant trades person and normal operation and use, in compliance with product instructions and information.

Warranty does not extend to product removal or re-installation. Warranty does not cover issues arising from incorrect installation.

---

# MOBILITY CARE MEMBERSHIP



## WARRANTY POLICY

Every Merits Mobility Scooter, purchased via a Goldfern Mobility reseller, now comes with a 3-year AA Mobility Care Membership.

Where a Nominated Mobility Scooter becomes immobile, the AA will attend, and, if unable to mobilise, either:

- Deliver the owner to their home address and arrange for the Mobility Scooter to be delivered to the owner's home address as soon as practicable.
- If the driver wishes to remain with the Mobility Scooter, the AA will arrange for recovery of both the driver and Mobility Scooter to the owner's home address, by means of a taxi van. This is free of charge to the driver as the cost of the taxi is charged back to the AA.

### TAXI ASSISTANCE - METROPOLITAN AREAS ONLY

Within metropolitan areas where Mobility Vans are available, the AA will make all the necessary arrangements to ensure both the driver and the scooter are transported to the owner's home address (or to any other address as specified by the driver at the time, providing the distance is no greater than to the owner's home address).

In order for the AA to provide service, the driver must remain with the Mobility Scooter until the AA Service Provider arrives.


The AA will not attend any Mobility Scooter that is immobile at the Owner's home address which requires repairs to be carried out.


---

## DISCLAIMER

All weights, dimensions, and measurements listed are approximate and may vary slightly. Minor variations may occur due to manufacturing tolerances, materials used, production processes, or design updates. These variations do not affect the overall performance or suitability of the product. Specifications are subject to change without notice.

**NOTES:**





06 357 4749



Orders@goldfernmobility.co.nz



208 John F Kennedy Drive, Milson, Palmerston North, NZ

**Goldfern**Mobility™

